

SiMPL Electronic Batch Record — Process Instruction (PI) Sheet

Sample Virus in 2000L Single-Use Bioreactor (SUB) — U.S. Use | **NEW (XStep-based EBR)**

Material / Product	8012xxx / Sample Virus in 2000L SUB
Process Order	_____
Batch #	_____
Control Recipe	_____
Master Recipe	_____
Plant / Bldg-Rm	_____
USDA License	XXXX BP MFG US 2000L SUB
Effective Date	_____

How to read this: this is the NEW digital EBR assembled from reusable XStep building blocks. Each step shows the original paper batch-record section it replaces ("Original BR: Section ____"), so it can be compared directly against *MABR — 100/2000L SUB Manufacturing Directions*. Recurring blocks (e.g. Daily Observations, Sample Submission, pH Adjustment, LAF Sanitization, Vessel Transfer) are shown once and re-used at each occurrence in the live recipe.

Old ↔ New Crosswalk

EBR Phase	New XStep	Original MABR Section
Batch Identification	SMPL: Virus Seed Identification	I
	SMPL: Cell Identification	II
Phase 1 — Production Vessel Preparation	SMPL: Vessel & Single-Use Bag Setup	III.A–B
	SMPL: Equipment Sterilization Verification	III.A + Equip. Steril. Record
	SMPL: Probe Identification	III.C
	SMPL: pH Probe Calibration	III.A, IV.G
	SMPL: Assembly & Line Connections	III.D–H, V.I
	SMPL: Controller & Gas Flow Verification	IV.A
	SMPL: Fill Volume Calculation	IV.C
Phase 2 — Media Preparation	SMPL: Component Addition	IV.B–C, IV.F
	SMPL: Medium Addition	IV.B (media fill)
	SMPL: Controller Setpoint & Alarm Configuration	IV.D, V.J, VI.A/F
	SMPL: DO Probe Calibration	IV.E
	SMPL: LAF Sanitization	IV.F, VI.J
	SMPL: Antifoam Addition	IV.F
	SMPL: Cell Receipt & Transfer	V.A–C
Phase 3 — Plant / Infection	SMPL: Daily Observations	V.D/M, VI (table)
	SMPL: Virus Seed Inoculation	V.D
	SMPL: Weighted-Average Titer	V.E
	SMPL: MOI Calculation	V.F
	SMPL: Thaw & Add Virus Seed	V.G
	SMPL: Actual MOI	V.H
	SMPL: Transfer Line Clip-Off	V.I
	SMPL: Inoculation Finish Time & SUB Label	V.K

	SMPL: Waste Inactivation	V.L
Phase 4 — Culture & Harvest	SMPL: Sample Submission	VI.K, XI.E/H/L
	SMPL: Harvest (Unifuge Centrifugation)	VI.B–I
	SMPL: Harvest Controller Shutdown	VI.F
	SMPL: Vessel Transfer	VI.I, X, XI.C
	SMPL: Cold Storage	VI.L
Phase 5 — SUB Filter Decontamination & Integrity	SMPL: SUB Filter Decontamination	VII.A
	SMPL: Filter Integrity Test	VII.B–C
Phase 6 — Pre-Inactivation	SMPL: Pre-Inactivation / Thaw	VIII.A–B
Phase 7 — Cyclization & Inactivation	SMPL: BEI Cyclization	IX.A–E
	SMPL: BEI Cyclization Hold	IX.C–D
	SMPL: Viral Inactivation	X.A–I
	SMPL: BEI Waste Neutralization	X.J
	SMPL: pH Adjustment	IX.E, X.F, XI.G
Phase 8 — Neutralization, Formulation & Dispense	SMPL: Neutralization	XI.A
	SMPL: Neutralization Hold	XI.B
	SMPL: Formulation	XI.D, XI.F
	SMPL: Bulk Dispense & Disposition	XI.C, I–J
	SMPL: Label Control & Reconciliation	VI.K, XI.K
Phase 9 — Closeout	SMPL: SAP & Excel Transaction Record	XI.M
	SMPL: Comments	Comments

Batch Identification

SMPL: Virus Seed Identification

Original BR: Section I

Instructions

Record the Master Seed identification and Last Testing Date below, then record each seed batch used with its Working Seed identification, passage, volume and titer (up to three batches). Max passage of seed is X+6.

Master Seed Identification *

Last Testing Date *

Data Input

#	Seed Batch	Seed Identification	Passage	Volume (mL)	Titer	Performed By *
1						

+ Add Row

Witness By *

SMPL: Cell Identification

Original BR: Section II

Instructions

Record the Master Cell Stock and Production Stock Cell identification. Max passage of cells must not exceed X+30. Record cell bioreactor number and testing dates.

Master Cell Stock Batch # *

Last Testing Date *

Production Stock Cell ID *

Batch # *

Batch Date *

Cell Bioreactor # *

Performed By *

Check By *

Phase 1 — Production Vessel Preparation

SMPL: Vessel & Single-Use Bag Setup

Original BR: Section III.A–B

Instructions

Install the single-use bag in the 2000L SUB per SOP XX-XXX-X. NOTE: Close 2 of the 4 SUB exhaust filters to prevent pressure accumulation (SUB is not designed for any pressure load); verify vent filters are not kinked. Record vessel and bag details.

Vessel Number *

Vessel Type2000L SUB

Single Use Bag Lot # *

Expiration Date *

Performed By *

Check By *

SMPL: Equipment Sterilization Verification

Original BR: Section III.A + Equip. Steril. Record

Instructions

Attach probe/equipment sterilization sticker to this MD. Attaching the sterilization sticker verifies that all packed sterilized equipment has been checked for integrity and is acceptable for use. Record sterilization run and expiry per phase (Media Prep / Inoculation / Harvest / Inactivation / Dispense).

Data Input					
#	Process Phase	Equipment ID	Sterilization Cycle / Run #	Sterilization Exp. Date	Performed By *
1					

+ Add Row

Witness By *

SMPL: Probe Identification

Original BR: Section III.C

Instructions

Record the probe number and sleeve serial number for each probe installed: DO probes, TruFluor DO, pH probes and TruFluor pH.

Data Input				
#	Probe Type	Probe No.	Sleeve Serial No.	Performed By *
1				

+ Add Row

Witness By *

SMPL: pH Probe Calibration

Original BR: Section III.A, IV.G

Instructions

Calibrate a bench-top pH meter per SOP XX-XXX. Adjust the bioreactor pH probe offset so it matches the offline reading. Verify the pH controller is in **manual mode** (no need to monitor pH for this process). Record the final pH at plant.

pH Probe No. *

Offline pH *

Probe Offset *

Final pH at Plant *

► Record

Date *

Time *

Performed By *

Check By *

SMPL: Assembly & Line Connections

Original BR: Section III.D–H, V.I

Instructions

Complete the following bioreactor connections per SOP XX-XXX-X, then sign off:

- Clamp off all lines to the bag except exhaust, overlay, the hollow spargers and probe lines.
- Connect overlay (Air Flow-1) and pressure sensor to the controller.
- Connect sparger 1 from the controller to the hollow sparger at the bottom of the bag.
- Insert DO and pH probes into the bag and connect their respective cables.
- Insert the temperature probe into the bag.

Performed By *

Check By *

Phase 2 — Media Preparation

SMPL: Controller & Gas Flow Verification

Original BR: Section IV.A

Instructions

From the controller's main screen, verify the following per SOP XX-XXX-X, then sign off:

- Air Flow-1 is flowing to the overlay.
- Air Flow-2, N2, CO2, O2 Flow-1 and O2 Flow-2 are flowing to the hollow sparger.
- Visually inspect the bottom of the bag for wrinkles.

Performed By *

Check By *

SMPL: Fill Volume Calculation

Original BR: Section IV.C

Instructions

Calculate the volume to fill the reactor before seed addition: $\text{Planned Final Volume} - \text{Seed Volume} = \text{Fill Volume}$. Record the calculation.

Standard Final Volume (L)

Planned Final Volume (kg) *

Seed Volume (kg) *

Fill Volume (kg)

Performed By *

Check By *

SMPL: Component Addition

Original BR: Section IV.B–C, IV.F

Instructions

Add components per the Inoculation Component Table. The component, manufacturer and standard planned amount are pre-filled; record the batch number, expiry date and amount added, then perform the SAP goods issue. Amounts are calculated per planned final volume.

Data Input

#	Component	Manufacturer	Amount	Batch Number	Exp. Date	Amount Added (L)	SAP	Performed By *
1	Excell Media	In House	QS to 2000L					

+ Add Row

Witness By *

SMPL: Medium Addition

Original BR: Section IV.B (media fill)

Instructions

Transfer 1100+/-55 kg of medium (PN: 8006775) into the bioreactor per SOP-0080692 or SOP-0080241. Specification for transfer amount is (1045 – 1155 kg). Post Addition Expiry Date is 5 days after medium charge (auto-calculated).

Data Input

#	Start	Transfer Start Date	Transfer Start Time	End	Transfer End Date	Transfer End Time	Amount Added	UoM	Post Addition Exp. Date	Performed By *
1	► Start	🕒	🕒	► End	🕒	🕒		kg		

+ Add Row

Witness By *

SMPL: Controller Setpoint & Alarm Configuration

Original BR: Section IV.D, V.J, VI.A/F

Instructions

Turn on the TCU connected to the SUB and set the controller setpoints per SOP XX-XXX-X. The standard set points are pre-filled below; record the actual value read back from the controller. Load the appropriate alarm file (Batch#_Alarms_ON at inoculation / Vessel_Process_Alarms_Off before harvest) as required.

Data Input

#	Parameter	Set Point	Actual	Performed By *
1	Temperature	2x°C		

+ Add Row

Witness By *

SMPL: DO Probe Calibration

Original BR: Section IV.E

Instructions

Once the vessel has reached the temperature set point, set the 0% and 100% calibration points for both electrochemical DO probes per SOP XX-XXX-X. Then set DO control to xx% on Auto, reduce stirrer to xx RPM, and zero the O2 totalizer. Enter the calibration values and stamp the completion time.

Data Input

#	DO Probe No.	Zero (nA)	Slope / Span (nA)		Date	Time	Performed By *
1				► Record	🕒	🕒	

+ Add Row

Witness By *

SMPL: LAF Sanitization

Original BR: Section IV.F, VI.J

Instructions

Sanitize the LAF with the standard disinfectant. Record the room, LAF, disinfectant batch and expiry date.

Disinfectant

Room #

LAF #

Disinfectant Batch #

Exp. Date

Performed By

SMPL: Antifoam Addition

Original BR: Section IV.F

Instructions

Aseptically, in a sterile vessel, add the appropriate amount of xxx mM L-Glutamine, Gentamycin, and antifoam according to the component table.

Performed By

Check By

Phase 3 — Plant / Infection

SMPL: Cell Receipt & Transfer

Original BR: Section V.A–C

Instructions

Receive cells from the cell prep lab and attach the cell prep information to this MD. Ensure the media is at 2x°C before transferring cells (up to 10 lb air pressure may be applied). Stamp the receipt and transfer times; record cell tank and sterilization expiry.

Data Input

#		Date Received	Time Received	Cell Tank #	Sterilization Exp. Date		Transfer Time	Performed By *
1	► Record	🕒	🕒			► Record	🕒	

+ Add Row

Witness By *

SMPL: Daily Observations

Original BR: Section V.D/M, VI (table)

Instructions

Sample the reactor daily and record bioreactor conditions in the Daily Observations table; stamp the sample date/time each day. Ensure parameters are in range (Temp 25–29°C, DO 20–100%, RPM 55). If out of range, contact the Supervisor and comment on any corrective action taken.

Data Input

#		Date	Time	Temp (°C)	RPM	pH / Probe	External pH	DO% / Probe	O2 Consumption (L)	CPE %	Cell Count Live	Cell Count Dead	Performed By *
1	► Record	🕒	🕒										

+ Add Row

Witness By *

SMPL: Virus Seed Inoculation

Original BR: Section V.D

Instructions

When the SUB reaches a cell density of approximately 1.0×10⁶ cells/mL, aseptically add the thawed virus seed to the bioreactor. Record the cell count and stamp the start and finish times of the addition.

Data Input

#	Cell Count (×10 ⁶ cells/mL)	Start	Start Time	End	End Time	Performed By *
1		► Start	🕒	► End	🕒	

+ Add Row

Witness By *

SMPL: Weighted-Average Titer

Original BR: Section V.E

Instructions

If multiple seed batches are used, calculate the weighted-average titer (refer to the Virus Seed Identification step). For each batch, record the Inv. Log of Titer and Seed Volume; the **Weighted Log** = Inv. Log of Titer × Seed Volume. The Final Inv. Log of Titer = (Total Weighted Log) ÷ (Total Volume) and is carried into the MOI Calculation.

Data Input

#	Seed Batch #	Inv. Log of Titer	Seed Volume (mL)	Weighted Log	Performed By *
1					

+ Add Row

Witness By *

SMPL: MOI Calculation

Original BR: Section V.F

Instructions

Calculate the amount of virus seed needed using the planned conditions for infection: (Reactor Fill Volume × Cell Density × MOI) ÷ Inv. Log Titer = Amount Seed Needed. Use the Final Inv. Log of Titer from the Weighted-Average Titer step if multiple seed batches were used.

Reactor Fill Volume (mL) *

Cell Density (cells/mL)

MOI

Inv. Log Titer *

Amount Seed Needed (mL)

Performed By *

Check By *

SMPL: Thaw & Add Virus Seed

Original BR: Section V.G

Instructions

Thaw the required amount of virus seed in a bath of disinfectant (the bath should be no warmer than 2x°C). Add the thawed seed to the bioreactor.

Performed By *

Check By *

SMPL: Actual MOI

Original BR: Section V.H

Instructions

After inoculation, calculate the Actual Plant MOI using the actual amount of seed added: (Actual Amount Seed Used × Inv. Log Titer) ÷ (Cell Density × Reactor Fill Volume) = Actual Plant MOI. Cell Density is the cell count at time of infection (carried from the Daily Observations / infection sample).

Actual Amount Seed Used (mL) *

Inv. Log Titer *

Cell Density (cells/mL) *

Reactor Fill Volume (mL) *

Actual Plant MOI

Performed By *

Check By *

SMPL: Transfer Line Clip-Off

Original BR: Section V.I

Instructions

Clipster off the lines used to transfer components, cells, and seed to the reactor, then sign off.

Performed By *

Check By *

SMPL: Inoculation Finish Time & SUB Label

Original BR: Section V.K

Instructions

Stamp the inoculation finish time and label the SUB as: Virus code / serial # / date.

Inoculation Finish Time *

► Record

SUB Label (Virus code / serial # / date) *

Performed By *

Check By *

SMPL: Waste Inactivation

Original BR: Section V.L

Instructions

Kill any remaining cells and seed with Hypochlor (per VSM 800.56) at 100 mL/L and discard. A minimum of 10 minutes disinfectant contact time is required before discard. Stamp the start/end times and record volumes and disinfectant.

Data Input

#	Start	Start Time	Vol. Seed Discarded (mL)	Vol. Cells Discarded (mL)	End	End Time	Disinfectant Batch #	Exp. Date	Performed By *
1	▶ Start	🕒			▶ End	🕒			

+ Add Row

Witness By *

Phase 4 — Culture & Harvest

SMPL: Sample Submission

Original BR: Section VI.K, XI.E/H/L

Instructions

Pull samples during processing and submit to QC. Record sample type, quantity/volume, storage condition and stamp the submission date. E.g. Retention 4×3mL freeze (-70°C); Bacteria 1×3mL store (4°C); Sterility / Densitometry as required.

Data Input

#	Sample Type	Qty × Volume	Storage Condition		Date Submitted	Performed By *
1				► Record	🕒	

+ Add Row

Witness By *

SMPL: Harvest (Unifuge Centrifugation)

Original BR: Section VI.B–I

Instructions

Connect the Unifuge to the harvest port. Harvest the antigen by passing 100% of the bioreactor volume through the Unifuge at a target 3 L/min, 3000 G. Record PBS used (SAP write-off), harvest tank weight, and stamp the start/finish times.

Data Input

#	Start	Start Time	Unifuge #	Module Lot / S/N	Harvest Tank #	Harvest Tank Steril. Exp.	SUB Vol. at Start (L)	PBS Batch #	PBS Exp. Date	PBS Amount (L)	SAP	End	End Time	Harvest Tank Weight (kg)	Performed By *
1	► Start	🕒										► End	🕒		

+ Add Row

Witness By *

SMPL: Harvest Controller Shutdown

Original BR: Section VI.F

Instructions

After harvesting approximately 60% of the SUB volume (~1200L for a 2000L batch), turn off the vent heaters, agitator, DO, and temperature control, then sign off.

Performed By *

Check By *

SMPL: Vessel Transfer

Original BR: Section VI.I, X, XI.C

Instructions

Aseptically transfer fluid between vessels (e.g., carboys to the inactivation bag/tank, or neutralized antigen to the dispense tank). Verify the scale reads zero before transfer where applicable. Record source, destination, volume transferred and stamp the start/end times.

Data Input

#	Start	Start Time	Source Vessel	Destination Vessel	Dest. Steril. Exp.	Volume Transferred (mL)	End	End Time	Performed By *
1	► Start	🕒					► End	🕒	

+ Add Row

Witness By *

SMPL: Cold Storage

Original BR: Section VI.L

Instructions

Seal carboys/pools with white tape and transfer to the freezer per SOP XX-XXX-X. First at ≤ -20°C then at ≤ -40°C. Stamp the date/time and record freezer and temperature.

Data Input

#	Container	Freezer #	Temp (°C)		Date	Time	Performed By *
1				► Record	🕒	🕒	

+ Add Row

Witness By *

Phase 5 — SUB Filter Decontamination & Integrity

SMPL: SUB Filter Decontamination

Original BR: Section VII.A

Instructions

Decontaminate the SUB filters by autoclave. Stamp the date and record the autoclave number and cycle number.

► Record

Date *

Autoclave # *

Cycle # *

Performed By *

SMPL: Filter Integrity Test

Original BR: Section VII.B–C

Instructions

Perform filter integrity testing on the SUB filters per SOP XX-XXX-X using the Millipore Integritest 5 instrument. Press Get Filter to populate the Filter ID and Filter Description. Attach integrity test results for each filter tested and record the result.

Data Input

#		Filter ID	Filter Description	Test Type	Result (Pass/Fail)		Date	Performed By *
1	► Get Filter					► Record	🕒	

+ Add Row

Witness By *

Phase 6 — Pre-Inactivation

SMPL: Pre-Inactivation / Thaw

Original BR: Section VIII.A–B

Instructions

Pull the harvest carboys from the freezer with ample time to thaw and warm to 3x°C. Pre-warm Sterile WFI for BEA/BEI cyclization to 3x°C the night before inactivation. Stamp the thaw/warm-up date/time; the target temperature is pre-filled.

Data Input

#	Item (Carboys / WFI)		Date	Time	Target Temp (°C)	Performed By *
1	Harvest Carboys	► Record	🕒	🕒	3x°C	

+ Add Row

Witness By *

Phase 7 — Cyclization & Inactivation

SMPL: BEI Cyclization

Original BR: Section IX.A–E

Instructions

Calculate volume of BEI needed: $\text{Antigen Volume} \div 9 = \text{Volume of BEI}$. The component and standard amount/mL are pre-filled; measure each component, combine BEA powder with WFI then add NaOH solution and adjust to final formulation volume with WFI. (Cyclization stir time and temperature are recorded in the BEI Cyclization Hold step; post-cyclization pH in the pH Adjustment step.) CAUTION: BEI is a potential mutagen — work in a ventilation hood / PAPR.

Data Input

#	Component	Manufacturer / Batch #	Exp. Date	Amount / mL	Volume to Add (mL)	SAP	Performed By *
1	BEA			0.0205 g			

+ Add Row

Witness By *

SMPL: BEI Cyclization Hold

Original BR: Section IX.C–D

Instructions

Stir the solution at $3\text{x}^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 1 hour to complete cyclization. Stamp the cyclization start and stop times and record the temperature at each. Use the BEI solution within x hours of the stop time.

Data Input

#	Start	Start Date	Start Time	Start Temp (°C)	End	End Time	Stop Temp (°C)	Performed By *
1	► Start	🕒	🕒		► End	🕒		

+ Add Row

Witness By *

SMPL: Viral Inactivation

Original BR: Section X.A–I

Instructions

Install the SUM bag in the reservoir tank; transfer the pre-warmed virus fluid and record the volume (scale value). Set Temperature $\text{XX}^{\circ}\text{C}$, Agitator XXX rpm . Add BEI solution through a $0.\text{xx} \mu\text{m}$ filter while mixing. Transfer to a second SUM bag, set parameters and enable alarms. Mix slowly at $3\text{x}^{\circ}\text{C} \pm 2^{\circ}\text{C}$ for $\sim\text{xx hr}$; stamp the inactivation start/stop times.

Data Input

#	Start	Start Time	Start Temp (°C)	SUM / Tank #	Bag Serial #	Steril. Exp. Date	Vol. Virus Fluid (mL)	Vol. BEI Added (mL)	End	Stop Time	Stop Temp (°C)	New Antigen Vol. (mL)	Performed By *
1	► Start	🕒							► End	🕒			

+ Add Row

Witness By *

SMPL: BEI Waste Neutralization

Original BR: Section X.J

Instructions

Neutralize the remaining BEI using 3M Na-thiosulfate: run it into the first SUM bag/tank through the drain tube, then place the bag in a yellow barrel and label it as Hazardous waste (if applicable). Record the BEI reconciliation and disinfectant details.

3M Na-Thiosulfate Batch # *

Exp. Date *

Volume of BEI Formulated (mL) *

Volume of BEI Added (mL) *

Volume of BEI Discarded (mL) *

Vol. 3M Na-Thio Added (mL) *

Performed By *

SMPL: pH Adjustment

Original BR: Section IX.E, X.F, XI.G

Instructions

Record the pH of the virus fluids and adjust to target (7.x ± 0.3) using HCl or NaOH as necessary. Record acid/base batch, expiry and amount added.

Initial pH *

Final pH *

HCl Batch # *

HCl Exp. Date *

NaOH Batch # *

NaOH Exp. Date *

Amount Added (mL) *

Performed By *

Check By *

Phase 8 — Neutralization, Formulation & Dispense

SMPL: Neutralization

Original BR: Section XI.A

Instructions

At the end of inactivation, add 1 mL of 3M Sodium Thiosulfate per XXX mL of inactivated virus fluids to yield a final concentration of 30 mM: Vol. of 3M Na-thio = Antigen Volume ÷ 100; New Antigen Volume = Antigen Volume + Vol. 3M Na-thio. Record the component in the Neutralization Component Table.

Data Input

#	Component	Batch Number	Exp. Date	Antigen Volume (mL)	Vol. 3M Na-Thio (mL)	New Antigen Vol. (mL)	SAP	Performed By *
1								

+ Add Row

Witness By *

SMPL: Neutralization Hold

Original BR: Section XI.B

Instructions

Allow the antigen to neutralize for 1 to 4 hours, mixing slowly at XXX rpm. Stamp the start and stop times and record the temperature at each.

Data Input

#	Start	Start Time	Start Temp (°C)	End	End Time	Stop Temp (°C)	Performed By *
1	► Start	🕒		► End	🕒		

+ Add Row

Witness By *

SMPL: Formulation

Original BR: Section XI.D, XI.F

Instructions

Add Gentamycin (final 0.000X mL/mL) and 10% Thimerosal (final 1:10,000). The component and standard final concentration are pre-filled. Pull pre- and post-Thimerosal samples and record each component addition and the resulting antigen volume. NOTE: samples must specify whether they contain Thimerosal.

Data Input

#	Component	Batch #	Exp. Date	Final Concentration	Volume Added (mL)	New Antigen Vol. (mL)	SAP	Performed By *
1	Gentamycin			0.000X mL/mL				

+ Add Row

Witness By *

SMPL: Bulk Dispense & Disposition

Original BR: Section XI.C, I–J

Instructions

Dispense the antigen into pools and store at ≤ -20°C. Label each pool with material #, batch #, pool #, pool volume and dispensing date. Record the bulk disposition and bottles used.

Data Input

#	Pool Number	Volume / Pool (mL)	Total Volume (mL)	Freezer #	Temp (°C)	Bottles Batch #	Bottles Used (ea)	Performed By *
1								

+ Add Row

Witness By *

SMPL: Label Control & Reconciliation

Original BR: Section VI.K, XI.K

Instructions

Record labels applied to carboys, pools and bottles (material #, batch #, pool #, volume, date). Verify all pools for label accuracy before sending to storage, and reconcile bottle/label usage.

Data Input

#	Item Labeled	Material #	Batch #	Quantity	Performed By *
1					

+ Add Row

Verified By *

SMPL: SAP & Excel Transaction Record

Original BR: Section XI.M

Instructions

Record completion of the batch SAP and Excel transactions (goods issues, write-offs, dispense postings). Record the reference and performer, and stamp the date for each system.

Data Input

#	Transaction System	Transaction / Reference		Date	Performed By *
1	SAP		► Record	🕒	

+ Add Row

Witness By *

SMPL: Comments

Original BR: Section Comments

Instructions

Record any comments, deviations, or corrective actions taken during the process. When proceeding with comments, the supervisor/designee is responsible for contacting the appropriate personnel.

Data Input

#	Comment / Deviation	Corrective Action	Performed By *
1			

+ Add Row

Witness By *